


# Accuro™

Immersion and external disinfection guidelines.



The following are approved for use with RIVANNA Accuro systems:

ULTRASOUND GELS	CLEANERS/DISINFECTANTS	HIGH-LEVEL DISINFECTANTS*
Aquagel	Cavi Wipes	Cidex*
Aquasonic 100	Transeptic Spray	
Aquasonic Clear	Sono wipes	
Clear Image	P.I. Spray II	
Clear Wave	Protex Wipes	
Scan	Sani-cloth HB Wipes	
Ultraphonic	Sani-cloth Plus Wipes	

\*Not to be used for immersion disinfection, external disinfection only.