



Immersion and external disinfection guidelines.

Accuro®



The following are approved for use with RIVANNA Accuro systems:

ULTRASOUND GELS	CLEANERS/DISINFECTANTS	HIGH-LEVEL DISINFECTANTS
Aquagel	Cavi Wipes	Sterris VPro Max2
Aquasonic 100	Transeptic Spray	Sterris VPro S2
Aquasonic Clear	Sono Wipes	Sterrad NX
Clear Image	P.I. Spray II	Sterrad 100NX
Clear Wave	Protex Wipes	Sterizone VP4
Scan	Sani-Cloth HB Wipes	
Ultraphonic	Sani-Cloth Plus Wipes	
	Sani-Cloth Bleach	

*Not to be used for immersion disinfection, external disinfection only.