



Immersion and external disinfection guidelines.

Accuro®



The following are approved for use with RIVANNA Accuro systems:

ULTRASOUND GELS	CLEANERS/DISINFECTANTS	HIGH-LEVEL DISINFECTANTS
Aquagel	Cavi Wipes	Cidex*
Aquasonic 100	Transeptic Spray	
Aquasonic Clear	Sono Wipes	
Clear Image	P.I. Spray II	
Clear Wave	Protex Wipes	
Scan	Sani-Cloth HB Wipes	
Ultraphonic	Sani-Cloth Plus Wipes	

*Not to be used for immersion disinfection, external disinfection only.